




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz still according to Anschütz-Thiele, straight,
with tempering jacket, flask capacity 50 ml,
sleeve/core NS 29/32, suction port GL 18, with
solid glass plug**

**€643.77
plus VAT &
Shipping**



Description

Intermediate receiver with straight thrust, with tempering jacket, graduated, with sleeve and core NS 29/32. Suction port with thread GL 18, NS stopcock with solid glass plug.
Made of DURAN® tube.

Additional Information

No.	LZ-1274
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3277
EAN	4058072310134
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32

