




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation receiver acc. to Anschütz-  
Thiele, straight, with tempering jacket, flask  
capacity 50 ml, sleeve/core NS 29/32, suction  
port GL 18, with PTFE plug**

**€653.99  
plus VAT &  
Shipping**



## Description

---

**Intermediate receiver** with straight thrust, with tempering jacket, graduated, with sleeve and core NS 29/32. Suction port with thread GL 18, NS stopcock with PTFE plug.  
Made of DURAN® tube.

## Additional Information

---

No.	LZ-1276
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3279
EAN	4058072310158
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32

