




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation receiver acc. to Anschütz-
Thiele, straight, with tempering jacket, flask
capacity 100 ml, sleeve/core NS 29/32, suction
port GL 18, with PTFE plug**

**€659.35
plus VAT &
Shipping**



Description

Intermediate receiver with straight thrust, with tempering jacket, graduated, with sleeve and core NS 29/32. Suction port with thread GL 18, NS stopcock with PTFE plug.
Made of DURAN® tube.

Additional Information

No.	LZ-1277
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3280
EAN	4058072310165
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32

