



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation receiver according to
Anschütz-Thiele, straight, flask capacity 100 ml,
sleeve/core NS 29/32, suction port GL 18, with
solid glass plug**

**€271.57
plus VAT &
Shipping**

Product Images



Description

Distilling masters according to Anschütz-Thiele; straight fore-stroke; graduated.

With NS sleeve and NS core, suction port with GL 18 thread, NS stopcock with solid glass plug.

Made of DURAN® tube.

Additional Information

No.	LZ-1264
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3267
EAN	4058072310035
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32

