



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation receiver acc. to Anschütz-
Thiele, slanted, flask capacity 50 ml,
sleeve/core NS 29/32, suction port GL 18, with
PTFE plug**

**€280.48
plus VAT &
Shipping**

Product Images



Description

Distilling trays according to Anschütz-Thiele; inclined forehearth; graduated.
With NS sleeve and NS core, suction port with GL 18 thread, NS stopcock with PTFE plug.
Made of DURAN® tube.

Additional Information

No.	LZ-1272
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3275
EAN	4058072310110
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	curved
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	29/32
Standard ground core	29/32

