



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation receiver acc. to Anschütz-  
Thiele, slanted, flask capacity 100 ml,  
sleeve/core NS 29/32, suction port GL 18, with  
PTFE plug**

**€287.51  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Distilling trays** acc. to Anschütz-Thiele; inclined forehearth; graduated.

With NS sleeve and NS core, suction port with GL 18 thread, NS stopcock with PTFE plug.

Made of DURAN® tube.

## Additional Information

---

No.	LZ-1273
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3276
EAN	4058072310127
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	curved
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	29/32
Standard ground core	29/32

