



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz still according to Bernauer, with 3 flasks 100 ml, sleeve top NS 19/26

**€113.42**  
plus VAT &  
Shipping

### Product Images



## Description

---

**Distilling masters** according to Bernauer (modified) with three cores NS 14/23 and 3 flasks.  
Suction nozzle with olive. Made of DURAN® tube.

## Additional Information

---

No.	LZ-1257
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3260
EAN	4058072309961
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Quantity pistons	3
Number of side necks	3
Alignment side necks	Oblique
Standard ground sleeve	19/26
Standard grinding for side necks Core	14/23

