



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distilling pattern according to Bernauer, sleeve top NS 14/23, bottom with 3 cores NS 14/23

€84.88
plus VAT &
Shipping

Product Images



Description

Distilling masters according to Bernauer (modified) with three cores NS 14/23; without flask.
Suction nozzle with olive. Made of DURAN® tube.

Additional Information

No.	LZ-1259
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3262
EAN	4058072309985
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	3
Standard ground sleeve	14/23
Alignment side necks	Oblique
Standard grinding for side necks Core	14/23

