



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz distilling pattern according to Bernauer, sleeve top NS 19/26, bottom with 3 cores NS 14/23

**€81.28**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Distilling masters** according to Bernauer (modified) with three cores NS 14/23; without flask.  
Suction nozzle with olive. Made of DURAN® tube.

## Additional Information

---

No.	LZ-1260
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3263
EAN	4058072309992
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	3
Standard ground sleeve	19/26
Alignment side necks	Oblique
Standard grinding for side necks Core	14/23

