



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz distilling pattern according to Bernauer, sleeve top NS 29/32, bottom with 3 cores NS 14/23

**€86.81**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Distilling masters** according to Bernauer (modified) with three cores NS 14/23; without flask.  
Suction nozzle with olive. Made of DURAN® tube.

## Additional Information

---

|                                       |                            |
|---------------------------------------|----------------------------|
| No.                                   | LZ-1261                    |
| Manufacturer (Brand)                  | Lenz                       |
| old neoLab article no.                | B-3264                     |
| EAN                                   | 4058072310004              |
| Transport temperature                 | Room temperature           |
| autoclavable at 121°C                 | Yes                        |
| Color                                 | transparent (clear)        |
| Material                              | Borosilicate glass         |
| Cooler features                       | with standard ground joint |
| Material                              | Borosilicate glass         |
| Number of side necks                  | 3                          |
| Standard ground sleeve                | 29/32                      |
| Alignment side necks                  | Oblique                    |
| Standard grinding for side necks Core | 14/23                      |

