



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation receiver n. Bredt, curved, sleeve top NS 14/13, middle sleeve-core connection NS 19/26, suction port GL 18, with 4 flasks 50 ml

**€203.61
plus VAT &
Shipping**

Product Images



Description

Distilling trays acc. to Bredt, bent

With four flasks NS 14/23. Suction port: GL 18.

Inclined thrust (105°)

Made of DURAN® tube.

Additional Information

No.	LZ-1245
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3248
EAN	4058072309848
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	bent 90°
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Quantity pistons	4
Number of side necks	4
Alignment side necks	Oblique 105°
Standard ground sleeve	19/26
Standard ground core	19/26
Standard grinding for side necks Core	14/23

