



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation receiver n. Bredt, curved, sleeve top NS 29/32, middle sleeve-core connection NS 29/32, suction port GL 18, with 4 flasks 100 ml

**€205.51
plus VAT &
Shipping**

Product Images



Description

Distilling trays according to Bredt, curved.

With four flasks NS 14/23. Suction port: GL 18.

Inclined thrust (105°)

Made of DURAN® tube.

Additional Information

No.	LZ-1246
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3249
EAN	4058072309855
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	bent 90°
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	4
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Quantity pistons	4
Standard ground sleeve	29/32
Standard ground core	29/32
Alignment side necks	Oblique 105°
Standard grinding for side necks Core	14/23

