



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation receiver n. Bredt, straight,  
sleeve top NS 14/23, middle sleeve-core  
connection NS 19/26, suction port GL 18, with 4  
flasks 50 ml**

**€203.61  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Distilling trays acc. to Bredt**, straight.  
With four flasks NS 14/23. Suction port: GL 18.  
Made of DURAN® tube.

## Additional Information

---

No.	LZ-1247
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3250
EAN	4058072309862
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	4
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Quantity pistons	4
Standard ground sleeve	19/26
Standard ground core	19/26
Alignment side necks	Oblique 105°
Standard grinding for side necks Core	14/23

