



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation forend bent, 105 °, sleeve/core NS 14/23, with drip tip

€18.52
plus VAT &
Shipping

Product Images



Description

Distillery foresets, bent (105°)
With NS sleeve and NS core, with drip tip.
Made of DURAN tube

Additional Information

No.	LZ-1217
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3220
EAN	4058072309565
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	14/23

