



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz distillation forend bent, 105 °, sleeve/core NS 29/32, with drip tip

**€25.16**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Distillery foresets, bent (105°)  
With NS sleeve and NS core, with drip tip.  
Made of DURAN tube

## Additional Information

---

No.	LZ-1220
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3223
EAN	4058072309596
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	29/32
Standard ground core	29/32

