



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation forend, bent 105 °, sleeve/core NS 14/23, with ventilation

€23.09
plus VAT &
Shipping

Product Images



Description

Distilling forcer, bent (105°), for high-boiling liquids.

With vent, with NS core, with drip tip.

Made of DURAN® tube.

Additional Information

No.	LZ-1225
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3228
EAN	4058072309640
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	14/23

