



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz distillation forend, bent 105 °, sleeve/core NS 29/32, with ventilation

**€29.38**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Distilling forcer**, bent (105°), for high-boiling liquids.

With vent, with NS core, with drip tip.

Made of DURAN® tube.

## Additional Information

---

No.	LZ-1228
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3231
EAN	4058072309671
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	29/32
Standard ground core	29/32

