



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation forend, straight, sleeve/core NS 14/23, with suction tube

€24.16
plus VAT &
Shipping

Product Images



Description

Distillery thrusts, straight. With NS sleeve and NS core.
With suction tube. Made of DURAN® tube.

Additional Information

| | |
|---|----------------------------|
| No. | LZ-1221 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3224 |
| EAN | 4058072309602 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Form | straight |
| Cooler features | with standard ground joint |
| Material | Borosilicate glass |
| Connection General / Tempering jacket / Cooler / Vacuum | Olive |
| Standard ground sleeve | 14/23 |
| Standard ground core | 14/23 |

