



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz distillation forend, straight, sleeve/core NS 24/29, with suction tube

**€24.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Distillery thrusts**, straight. With NS sleeve and NS core.  
With suction tube. Made of DURAN® tube.

## Additional Information

---

No.	LZ-1223
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3226
EAN	4058072309626
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Standard ground sleeve	24/29
Standard ground core	24/29

