



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation forend, straight, sleeve/core NS 29/32, with suction tube

€26.06
plus VAT &
Shipping

Product Images



Description

Distillery thrusts, straight. With NS sleeve and NS core.

With suction tube. Made of DURAN® tube.

Additional Information

No.	LZ-1224
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3227
EAN	4058072309633
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Standard ground sleeve	29/32
Standard ground core	29/32

