



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation forend with long outlet, bent, sleeve NS 14/23, tube length 190 mm

€19.72
plus VAT &
Shipping

Product Images



Description

Distillery bumpers with long outlet, bent (105°); with NS sleeve.

Made of DURAN® tube.

Additional Information

No.	LZ-1213
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3216
EAN	4058072309527
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Length	190 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23

