



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation forend with long drain, bent,  
sleeve NS 24/29, tube length 200 mm**

**€21.56  
plus VAT &  
Shipping**

## Product Images



## Description

**Distillery bumpers** with long outlet, bent (105°); with NS sleeve.  
Made of DURAN® tube.

## Additional Information

|                        |                            |
|------------------------|----------------------------|
| No.                    | LZ-1215                    |
| Manufacturer (Brand)   | Lenz                       |
| old neoLab article no. | B-3218                     |
| EAN                    | 4058072309541              |
| Transport temperature  | Room temperature           |
| autoclavable at 121°C  | Yes                        |
| Color                  | transparent (clear)        |
| Material               | Borosilicate glass         |
| Form                   | curved 105°                |
| Length                 | 200 mm                     |
| Cooler features        | with standard ground joint |
| Material               | Borosilicate glass         |
| Standard ground sleeve | 24/29                      |

