



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz vacuum advance bent, with threaded connection GL 18, sleeve/core NS 29/32

**€46.32**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Vacuum thrusts**, bent, with threaded connection

According to DIN 12594, angle (105°), with NS sleeve and NS core.

With drip tip; vacuum connection: olive

## Additional Information

---

No.	LZ-1244
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3247
EAN	4058072309831
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	29/32
Standard ground core	29/32
standardized according to DIN standard (DIN)	DIN 12594

