



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz vacuum pusher bent, with olive, sleeve/core NS 14/23

€32.49
plus VAT &
Shipping

Product Images



Description

Vacuum thrusts, bent, with olive

According to DIN 12594, angle (105°); with NS sleeve and NS core.

With drip tip; vacuum connection: olive

Additional Information

No.	LZ-1237
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3240
EAN	4058072309763
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	14/23
standardized according to DIN standard (DIN)	DIN 12594

