



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz vacuum advance bent, with olive, sleeve/core NS 29/23

**€39.17**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Vacuum thrusts**, bent, with olive

According to DIN 12594, angle (105°); with NS sleeve and NS core.

With drip tip; vacuum connection: olive

## Additional Information

---

No.	LZ-1240
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3243
EAN	4058072309794
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	29/32
Standard ground core	29/32
standardized according to DIN standard (DIN)	DIN 12594

