



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz vacuum advance straight, with threaded connection GL 18, sleeve/core NS 14/23

€34.81
plus VAT &
Shipping

Product Images



Description

Vacuum thrusts, straight, with threaded connection

With NS sleeve and NS core; with drip tip.

Vacuum connection: thread GL 18.

Additional Information

No.	LZ-1233
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3236
EAN	4058072309725
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	14/23
Standard ground core	14/23
standardized according to DIN standard (DIN)	DIN 12594

