



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz vacuum pusher straight, with olive, sleeve/core NS 24/29

**€29.46**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Vacuum plungers**, straight, with olive according to DIN 12594, with NS sleeve and NS core, with drip tip.  
Vacuum connection: olive

## Additional Information

---

No.	LZ-1231
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3234
EAN	4058072309701
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	24/29
Standard ground core	24/29
standardized according to DIN standard (DIN)	DIN 12594

