



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz foreshaft for distillation receiver acc. to
Bredt, bent, sleeve top NS 29/32, sleeve bottom
NS 29/32**

**€59.93
plus VAT &
Shipping**

Product Images



Description

Bendsfor Bredt distillation masters.

With two NS sleeves, suction connection: GL 18.

Curved forejoints: 105°C

Made of DURAN® tube.

Additional Information

No.	LZ-1251
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3254
EAN	4058072309909
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	curved 105°
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	29/32

