



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz foreshaft for distillation receiver acc. to  
Bredt, straight, sleeve top NS 14/23, sleeve  
bottom NS 19/26**

**€47.84  
plus VAT &  
Shipping**

**Product Images**



## Description

---

**Foresets**, straight for distilling masters according to Bredt.  
With two NS sleeves, suction connection: GL 18.  
Made of DURAN® tube.

## Additional Information

---

No.	LZ-1252
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3255
EAN	4058072309916
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	14/23 19/26

