



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz foreshaft for distillation receiver acc. to Bredt, straight, sleeve top NS 29/32, sleeve bottom NS 29/32

**€52.42**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Foresets**, straight for distillation masters according to Bredt.  
With two NS sleeves, suction connection: GL 18.  
Made of DURAN® tube.

## Additional Information

---

No.	LZ-1253
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3256
EAN	4058072309923
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Standard ground sleeve	29/32

