



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler for extractors, with olives,
core NS 34/35, extractor 70 ml, cooling surface
160 cm²**

**€104.69
plus VAT &
Shipping**



Description

Dimroth coolers for extractors with standard ground joint for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9-12 mm are required for this.

Additional Information

No.	LZ-1186
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3188
EAN	4058072309251
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	70 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	160 cm ²
Material olives	Borosilicate glass
Standard ground core	34/35
Type cooler	Dimroth cooler

