



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler for extractors, with olives,
core NS 60/46, extractor 300 – 500 ml, cooling
surface 680 cm²**

**€158.10
plus VAT &
Shipping**

Product Images



Description

Dimroth coolers for extractors with standard ground joint for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9-12 mm are required for this.

Additional Information

No.	LZ-1188
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3190
EAN	4058072309275
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	500 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	680 cm ²
Standard ground core	60/46
Type cooler	Dimroth cooler
Material olives	Borosilicate glass

