



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz Dimroth cooler for extractors, with olives,  
core NS 71/51, extractor 1,000 ml, cooling surface  
930 cm<sup>2</sup>**

**€169.28  
plus VAT &  
Shipping**



## Description

---

Dimroth coolers for extractors with standard ground joint for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9-12 mm are required for this.

## Additional Information

---

No.	LZ-1189
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3191
EAN	4058072309282
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	1000 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	930 cm <sup>2</sup>
Standard ground core	71/51
Type cooler	Dimroth cooler
Material olives	Borosilicate glass

