



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz Dimroth cooler for extractors, with GL connections, core NS 45/40, extractor 100 - 250 ml, water connection GL 14, cooling surface 370 cm²

**€124.51
plus VAT &
Shipping**

Product Images

Description

Dimroth coolers for extractors with standard ground joint for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are an indispensable part of the chemical synthesis laboratory nowadays. Design: Water connection with GL 14 thread (with two pre-mounted, unscrewable plastic olives)

Additional Information

No.	LZ-1192
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3194
EAN	4058072309312
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	250 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	370 cm ²
Standard ground core	45/40
Type cooler	Dimroth cooler
Material olives	Polypropylene (PP)

