



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler with olives, sleeve/core NS
19/26, shell length 160 mm, cooling surface 140
cm²**

**€89.06
plus VAT &
Shipping**

Product Images



Description

Dimroth standard ground joint condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9-12 mm are required for this.

Additional Information

No.	LZ-1168
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3170
EAN	4058072309077
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 160
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	140 cm ²
Material olives	Borosilicate glass
Standard ground sleeve	19/26
Standard ground core	19/26
Type cooler	Dimroth cooler

