



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler with GL connections,
sleeve/core NS 14/23, shell length 160 mm,
cooling surface 140 cm²**

**€104.12
plus VAT &
Shipping**

Product Images



Description

Dimroth coolers with standard ground joint for the condensation of gases and vapors in chemical synthesis.

These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory.

Design:

Water connection with GL 14 thread (with two pre-mounted, unscrewable plastic olives)

Additional Information

No.	LZ-1169
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3171
EAN	4058072309084
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 160
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	140 cm ²
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Dimroth cooler
Material olives	Polypropylene (PP)

