



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler with olives, sleeve/core NS
29/32, jacket length 250 mm, cooling surface
440 cm²**

**€111.61
plus VAT &
Shipping**



Description

Dimroth standard ground joint condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9-12 mm are required for this.

Additional Information

No.	LZ-1174
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3176
EAN	4058072309138
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 250
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	440 cm ²
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Dimroth cooler
Material olives	Borosilicate glass
standardized according to DIN standard (DIN)	DIN 12591

