



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz Dimroth cooler with GL connections,
sleeve/core NS 45/40, shell length 400 mm,
cooling surface 1380 cm²**

**€203.13
plus VAT &
Shipping**

Product Images



Description

Dimroth standard ground joint condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Water connection with GL 14 thread (with two pre-mounted, unscrewable plastic olives)

Additional Information

No.	LZ-1184
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3186
EAN	4058072309237
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	1380 cm ²
Standard ground sleeve	45/40
Standard ground core	45/40
Type cooler	Dimroth cooler
Material olives	Polypropylene (PP)

