



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz pycnometer attachment with wide-neck steep-bottomed bottle, capacity 1,000 ml

€96.90
plus VAT &
Shipping

Product Images



Description

Pycnometer attachments with wide-necked steep-breasted bottle, with NS ground joint.
For determining the density (bulk density) of bituminous masses and mineral mixtures.

Additional Information

No.	LZ-1575
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3578
EAN	4058072313135
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass 3.3
Volume	1000 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	60/37
Standard ground core	60/37

