



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation attachment n. Claisen, 2 sleeves parallel NS 14/23, core vertical NS 24/29, core oblique NS 19/26

€51.86
plus VAT &
Shipping

Product Images



Description

Distillation tops according to Claisen (B-2934 to B-2939)

With two parallel sleeves NS 14/23 and two NS cores (drain angle: 75°)

For distillation thermometers with an installation length of 55 mm

Made of DURAN® tube

Additional Information

No.	LZ-0936
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2936
EAN	4058072306762
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	2
Standard ground sleeve	14/23
Standard ground core	24/29
Alignment side necks	Drain angle 75
Standard grinding for side necks Core	19/26

