



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation attachment n. Claisen, 2 sleeves parallel NS 14/23, core vertical NS 29/32, core oblique NS 14/23

**€44.36
plus VAT &
Shipping**

Product Images



Description

Distillation tops according to Claisen (B-2934 to B-2939)

With two parallel sleeves NS 14/23 and two NS cores (drain angle: 75°)

For distillation thermometers with an installation length of 55 mm

Made of DURAN® tube

Additional Information

No.	LZ-0938
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2938
EAN	4058072306786
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Number of side necks	2
Alignment side necks	Drain angle 75
Standard ground sleeve	14/23
Standard ground core	29/32
Standard grinding for side necks Core	14/23

