



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz dosing attachment for tipping machine, dosing capacity 3 ml

€44.17
plus VAT &
Shipping

Product Images



Description

Dosing attachments for automatic tilting machines with core 29/32.

Additional Information

No.	LZ-1456
Manufacturer (Brand)	Lenz
old neolab article no.	B-3459
EAN	4058072311957
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass 3.3
Volume	3 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	29/32
Standard ground core	29/32

