



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction attachment for specific light
extraction agents, extractor 500 ml, cooler
connection NS 60/46, with distribution tube**

**€315.46
plus VAT &
Shipping**

Product Images



Description

Extraction attachments for specific light extractants
For extraction from liquids with specific light extractants.
With insert filter (porosity: 0)

Additional Information

No.	LZ-1384
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3387
EAN	4058072311230
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	60/46
Standard ground core	60/46

