



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction attachment acc. to Soxhlet,
extractor 150 ml, cooler connection NS 45/40,
with NS glass tap**

**€142.07
plus VAT &
Shipping**

Product Images



Description

Extraction attachments in addition to the DIN standard.
With NS glass stopcock.

Additional Information

No.	LZ-1351
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3354
EAN	4058072310905
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	150 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	45/40
Standard ground core	45/40

