



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction attachment according to
Soxhlet, extractor 100 ml, cooler connection NS
45/40, without stopcock**

**€103.14
plus VAT &
Shipping**



Description

Extraction attachments according to Soxhlet

According to DIN 12602 and in addition to the DIN standard.

Without NS glass tap

Additional Information

No.	LZ-1342
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3345
EAN	4058072310813
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	45/40
Standard ground core	45/40

