



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction attachment according to
Soxhlet, extractor 200 ml, cooler connection NS
45/40, without stopcock**

**€109.29
plus VAT &
Shipping**



Description

Extraction attachments according to Soxhlet

In addition to the DIN standard.

Without NS glass stopcock

Additional Information

No.	LZ-1344
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3347
EAN	4058072310837
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	200 ml
Form	straight
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	45/40
Standard ground core	45/40

