



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz high performance cooler with 4 cooling coils, ball and socket KF 50, water connection balls KF 15

€1,336.10
plus VAT &
Shipping

Product Images



Description

High-performance cooler with four cooling coils
Cooling surface: 2,350 cm², shell length: 300 mm,
made of DURAN® tube

Additional Information

No.	LZ-1208
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3211
EAN	4058072309473
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 300
Material	Borosilicate glass
Connection cooling medium	KF 15
Spherical cut ball	KF 50
Cooling surface	3000 cm ²
Type cooler	High performance cooler with four cooling coils
Spherical cut bowl	KF 50
Material olives	Borosilicate glass DURAN®

