



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz intensive cooler with olives, sleeve/core NS
24/29, jacket length 250 mm, cooling surface
400 cm²**

**€135.33
plus VAT &
Shipping**

Product Images



Description

Intensive condenser (double-walled condenser) with standard ground joint for the condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives

Additional Information

| | |
|---------------------------|----------------------------|
| No. | LZ-1199 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3201 |
| EAN | 4058072309381 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Length | 0 250 |
| Cooler features | with standard ground joint |
| Material | Borosilicate glass |
| Connection cooling medium | Glass olive |
| Cooling surface | 400 cm ² |
| Standard ground sleeve | 24/29 |
| Standard ground core | 24/29 |
| Type cooler | Intensive cooler |
| Material olives | Borosilicate glass |

