



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz intensive cooler with olives, sleeve/core NS  
29/32, jacket length 400 mm, cooling surface  
650 cm<sup>2</sup>**

**€131.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

Intensive condenser (double-walled condenser) with standard ground joint for the condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives

## Additional Information

---

No.	LZ-1203
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3205
EAN	4058072309428
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	650 cm <sup>2</sup>
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Intensive cooler
Material olives	Borosilicate glass
standardized according to DIN standard (DIN)	DIN 12593

