



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz column n. Hempel with vacuum jacket,
filling height 200 mm, sleeve/core NS 14/23**

**€164.73
plus VAT &
Shipping**

Product Images



Description

Columns according to Hempel with fused-on vacuum jacket.
With NS sleeve and NS core. For filling with packing (e.g. glass beads).
Made of DURAN® tube.

Additional Information

No.	LZ-1061
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3061
EAN	4058072308018
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Height	200 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	14/23

